

Product Name: MCFD2 Mouse Monoclonal Antibody
Catalog #: AMM86057



Summary

Production Name	MCFD2 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	MCFD2
Alternative Names	Multiple coagulation factor deficiency protein 2, Neural stem cell-derived neuronal survival protein, MCFD2, SDNSF
Gene ID	90411.0
SwissProt ID	Q8NI22. This MCFD2 monoclonal antibody is generated from mouse immunized with MCFD2 recombinant protein.

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	16.4kDa

Background


Product Name: MCFD2 Mouse Monoclonal Antibody
Catalog #: AMM86057



The MCFD2-LMAN1 complex forms a specific cargo receptor for the ER-to-Golgi transport of selected proteins. Plays a role in the secretion of coagulation factors.

Research Area

Image Data



Anti-MCFD2 Antibody at 1:1000 dilution + 293T/17 whole cell lysate

Note

For research use only.