
Product Name: THEM2 Mouse Monoclonal Antibody**Catalog #: AMM86059**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	15.0kDa

Antigen Information


Gene Name	THEM2
Alternative Names	Acyl-coenzyme A thioesterase 13, Acyl-CoA thioesterase 13, 312-, Thioesterase superfamily member 2, Acyl-coenzyme A thioesterase 13, N-terminally processed, ACOT13, THEM2
Gene ID	55856.0
SwissProt ID	Q9NPJ3
Immunogen	This THEM2 monoclonal antibody is generated from mouse immunized with THEM2 recombinant protein.

Background

THEM2 belongs to the UPF0152 (paal) family. Its exact function is still undetermined.

Research Area

Image Data



All lanes : Anti-THEM2 Antibodyat 1:1000 dilution