
Product Name: CBFA2T3 Mouse Monoclonal Antibody**Catalog #: AMM86062**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	71.2kDa

Antigen Information


Gene Name	CBFA2T3 Protein CBFA2T3, MTG8-related protein 2, Myeloid translocation gene on chromosome 16
Alternative Names	protein, hMTG16, Zinc finger MYND domain-containing protein 4, CBFA2T3, MTG16, MTGR2, ZMYND4
Gene ID	863.0
SwissProt ID	O75081
Immunogen	This CBFA2T3 antibody is generated from a mouse immunized with a recombinant protein between 418-653 amino acids from human CBFA2T3.

Background

Functions as a transcriptional repressor. Regulates the proliferation and the differentiation of erythroid progenitors by repressing the expression of TAL1 target genes. Plays a role in granulocyte differentiation.

Research Area

Image Data



All lanes : Anti-CBFA2T3 Antibody at 1:500-1:2000 dilution