

Product Name: GOLGA5 Mouse Monoclonal Antibody**Catalog #: AMM86063**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	83.0kDa

Antigen Information

Gene Name	GOLGA5
Alternative Names	Golgin subfamily A member 5, Cell proliferation-inducing gene 31 protein, Golgin-84, Protein Ret-II, RET-fused gene 5 protein, GOLGA5, RETII, RFG5
Gene ID	9950.0
SwissProt ID	Q8TBA6
Immunogen	This GOLGA5 antibody is generated from a mouse immunized with a recombinant protein between 1-300 amino acids from human GOLGA5.


Background

Involved in maintaining Golgi structure. Stimulates the formation of Golgi stacks and ribbons. Involved in intra-Golgi

retrograde transport.

Research Area

Image Data



All lanes : Anti-GOLGA5 Antibody at 1:2000 dilution