
Product Name: WBP2 Mouse Monoclonal Antibody**Catalog #: AMM86077**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	28.1kDa


Antigen Information

Gene Name	WBP2
Alternative Names	WW domain-binding protein 2, WBP-2, WBP2
Gene ID	23558.0
SwissProt ID	Q969T9
Immunogen	This WBP2 antibody is generated from a mouse immunized with a recombinant protein between 1-261 amino acids from the C-terminal region of human WBP2.

Background

Research Area

Image Data



All lanes : Anti-WBP2 Antibody (C-term) at 1:2000 dilution