

Product Name: RAB23 Mouse Monoclonal Antibody
Catalog #: AMM86080



Summary

Production Name	RAB23 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,ICC,FC
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	RAB23
Alternative Names	Ras-related protein Rab-23, RAB23
Gene ID	51715.0
SwissProt ID	Q9ULC3.Purified His-tagged RAB23 protein(Fragment) was used to produced this monoclonal antibody.

Application

Dilution Ratio	WB 1:1000-1:2000,ICC 1:25-1:50,FC 1:25-1:50
Molecular Weight	26.7kDa

Background


Product Name: RAB23 Mouse Monoclonal Antibody
Catalog #: AMM86080



The protein encoded by this gene belongs to the small GTPase superfamily, Rab family. It may be involved in small GTPase mediated signal transduction and intracellular protein transportation. Alternative splicing occurs at this locus and two transcript variants encoding the same protein have been identified.

Research Area

Image Data



All lanes : Anti-RAB23 Antibody at 1:2000 dilution

Note

For research use only.