

---

**Product Name: VILIP1 Mouse Monoclonal Antibody****Catalog #: AMM86095**

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	WB
<b>Reactivity</b>	Human, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG1
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide.
<b>Purification</b>	Affinity Purification

**Application**

<b>Dilution Ratio</b>	WB 1:1000-1:2000
<b>Molecular Weight</b>	22.1kDa

**Antigen Information**


<b>Gene Name</b>	VILIP1
<b>Alternative Names</b>	Visinin-like protein 1, VILIP, VLP-1, Hippocalcin-like protein 3, HLP3, VSNL1, VISL1
<b>Gene ID</b>	7447.0
<b>SwissProt ID</b>	P62760
<b>Immunogen</b>	This VILIP1 antibody is generated from a mouse immunized with a recombinant protein of human VILIP1.

**Background**

Regulates (in vitro) the inhibition of rhodopsin phosphorylation in a calcium-dependent manner.

## Research Area

## Image Data



All lanes : Anti-VILIP1 Antibody at 1:4000 dilution