

Product Name: DRP-2 Mouse Monoclonal Antibody**Catalog #: AMM86096**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:2000-1:8000
Molecular Weight	62.3kDa


Antigen Information

Gene Name	DRP-2
Alternative Names	Dihydropyrimidinase-related protein 2, DRP-2, Collapsin response mediator protein 2, CRMP-2, N2A3, Unc-33-like phosphoprotein 2, ULIP-2, DPYSL2, CRMP2, ULIP2
Gene ID	1808.0
SwissProt ID	Q16555
Immunogen	This DRP-2 antibody is generated from mice immunized with a KLH conjugated synthetic peptide between 463-492 amino acids from human DRP-2.

Background

Research Area

Image Data



All lanes : Anti-DRP-2 Antibody at 1:8000 dilution