

---

**Product Name:** human CD25 Mouse Monoclonal Antibody**Catalog #:** AMM86105

For research use only.

**Summary**

<b>Description</b>	Mouse monoclonal Antibody
<b>Host</b>	Mouse
<b>Application</b>	FC
<b>Reactivity</b>	Human
<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Mouse IgG2a
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Concentration</b>	1mg/ml
<b>Storage</b>	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
<b>Shipping</b>	Ice bags
<b>Buffer</b>	Purified antibody in TBS with 0.05% sodium azide.
<b>Purification</b>	Affinity Purification

**Application**


<b>Dilution Ratio</b>	FC 1:50-1:100
<b>Molecular Weight</b>	/

**Antigen Information**

<b>Gene Name</b>	human CD25
<b>Alternative Names</b>	human CD25
<b>Gene ID</b>	
<b>SwissProt ID</b>	unknow
<b>Immunogen</b>	

**Background****Research Area**

## Image Data



Overlay histogram showing PHA-stimulated human peripheral blood lymphocytes stained with CD25 antibody (green line). The cells were incubated in 2% bovine serum albumin to block non-specific protein-protein interactions followed by the antibody (1:50 dilution) for 60min at 37°C. The secondary antibody used was Goat Anti-Mouse IgG, DyLight® 488 Conjugated Highly Cross-Adsorbed (OJ192088) at 1/200 dilution for 40min at 37°C. Isotype control antibody (blue line) was mouse IgG2a (1 µg/1x10<sup>6</sup> cells) used under the same conditions. Acquisition of >10,000 events was performed.