
Product Name: PRDM12 Mouse Monoclonal Antibody**Catalog #: AMM86111**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in TBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:500-1:2000
Molecular Weight	40.4kDa


Antigen Information

Gene Name	PRDM12
Alternative Names	PR domain zinc finger protein 12, 211-, PR domain-containing protein 12, PRDM12, PFM9
Gene ID	59335.0
SwissProt ID	Q9H4Q4
Immunogen	This PRDM12 antibody was raised using purified recombinant GST fusion protein encoding human PRDM12.

Background

Research Area

Image Data



Anti-PRDM12 Antibody at 1:1000 dilution + 293 whole cell lysate