

Summary

Production Name	ITGA11 Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	ITGA11
Alternative Names	Integrin alpha-11, ITGA11
Gene ID	22801.0
SwissProt ID	Q9UKX5. This ITGA11 antibody is generated from a mouse immunized with a KLH conjugated synthetic peptide between 675-880 amino acids from human ITGA11.

Application

Dilution Ratio	WB:1:8000
Molecular Weight	133.5kDa

Background

Integrin alpha-11/beta-1 is a receptor for collagen.


Product Name: ITGA11 Mouse Monoclonal Antibody
Catalog #: AMM86119



Research Area

PI3K-Akt signaling pathway

Image Data



All lanes : Anti-ITGA11 Antibody at 1:8000 dilution

Note

For research use only.