

Product Name: GSTA1 Mouse Monoclonal Antibody**Catalog #: AMM86122**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:2000-1:5000
Molecular Weight	25.6kDa

Antigen Information


Gene Name	GSTA1
Alternative Names	Glutathione S-transferase A1, 2.5.1.18, GST HA subunit 1, GST class-alpha member 1, GST-epsilon, GSTA1-1, GTH1, Glutathione S-transferase A1, N-terminally processed, GSTA1
Gene ID	2938.0
SwissProt ID	P08263
Immunogen	This GSTA1 antibody is generated from a mouse immunized with a recombinant protein from the human region of human GSTA1.

Background

Conjugation of reduced glutathione to a wide number of exogenous and endogenous hydrophobic electrophiles.

Research Area

Image Data



All lanes : Anti-GSTA1 Antibody at 1:5000 dilution