

Product Name: SET07 Mouse Monoclonal Antibody
Catalog #: AMM86128



Summary

Production Name	SET07 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ICC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide.
Purification	Affinity Purification

Immunogen

Gene Name	SET07 N-lysine methyltransferase SETD8, 211-, H4-K20-HMTase SETD8, Histone-lysine N-methyltransferase SETD8, Lysine N-methyltransferase 5A, PR/SET domain-containing protein 07, PR-Set7, PR/SET07, SET domain-containing protein 8, SETD8, KMT5A, PRSET7, SET07, SET8
Alternative Names	
Gene ID	387893.0
SwissProt ID	Q9NQR1.This SET07 antibody was raised using purified recombinant GST fusion protein encoding the N-terminal region of human SET07.

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:100-1:500,ICC 1:20-1:50
Molecular Weight	42.9kDa


Product Name: SET07 Mouse Monoclonal Antibody
Catalog #: AMM86128



Background

Research Area

Image Data



All lanes :SET07 Antibody at 1:2000 dilution

Note

For research use only.