

**Product Name: Phosphotyrosine(2C1) Mouse
Monoclonal Antibody
Catalog #: AMM86139**



Summary

Production Name	Phosphotyrosine(2C1) Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Other

Performance

Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% BSA and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Phosphotyrosine(2C1)
Alternative Names	Phosphotyrosine
Gene ID	
SwissProt ID	/Recombinant Protein

Application

Dilution Ratio	WB:1:1000-1:3000 ,IHC:1:50-1:200
Molecular Weight	/


**Product Name: Phosphotyrosine(2C1) Mouse
Monoclonal Antibody
Catalog #: AMM86139**




Background

Research Area

Image Data



Western blot analysis of 1) Jurkat, 2) 293T, 3) Rat Brain Tissue, 4) Mouse Brain Tissue with Phosphotyrosine Mouse mAb diluted at 1:2,000.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Phosphotyrosine Mouse mAb diluted at 1:200.

Note

For research use only.