

Product Name: CaMKII β / γ / δ (Phospho Thr287)(4H2) Mouse Monoclonal Antibody
Catalog #: AMM86140

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Phosphorylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:100-1:200
Molecular Weight	50kDa


Antigen Information

Gene Name	CaMKII β / γ / δ (Phospho Thr287)(4H2)
Alternative Names	CaMKII β / γ / δ (Phospho Thr287)(4H2)
Gene ID	816;
SwissProt ID	Q13554
Immunogen	Synthetic Peptide


Background

Research Area

Image Data



Immunohistochemical analysis of paraffin-embedded Rat Brain Tissue using CaMKIIβ/γ/δ(Phospho Thr287)(4H2) antibody diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Mouse Brain Tissue using CaMKIIβ/γ/δ(Phospho Thr287)(4H2) Mouse Monoclonal Antibody diluted at 1:200.