

Product Name: Cytokeratin 6 Mouse Monoclonal Antibody
Catalog #: AMM86153

Summary

Production Name	Cytokeratin 6 Mouse Monoclonal Antibody
Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB,IHC,ICC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Mouse IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Cytokeratin 6
Alternative Names	Cytokeratin 6 Mouse
Gene ID	
SwissProt ID	P02538/P04259/P48668.

Application

Dilution Ratio	WB 1:1000-1:2000,IHC 1:50-1:300,ICC 1:100-1:200
Molecular Weight	46-58kDa

Background


Product Name: Cytokeratin 6 Mouse Monoclonal Antibody
Catalog #: AMM86153



Cytokeratins are proteins of keratin-containing intermediate filaments found in the intracytoplasmic cytoskeleton of epithelial tissue.

Research Area

Image Data



Western blot analysis of 1) HeLa, 2) HepG2 , using AMM86153 diluted at 1:2,000.

Note

For research use only.