
Product Name: HAO1 Mouse Monoclonal Antibody**Catalog #: AMM86161**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:500-1:2000
Molecular Weight	41kDa

Antigen Information

Gene Name	HAO1
Alternative Names	HAO1
Gene ID	54363;
SwissProt ID	Q9UJM8
Immunogen	Recombinant Protein

Background**Research Area**

Image Data



Western blot analysis of 1) Mouse Liver Tissue, 2) Rat Liver Tissue using HAO1 Monoclonal Antibody.