
Product Name: Rabbit IgG (Fc specific, Heavy Chain) Mouse Monoclonal Antibody
Catalog #: AMM86162

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	WB 1:2000-1:5000
Molecular Weight	55kDa


Antigen Information

Gene Name	Rabbit IgG (Fc specific, Heavy Chain)
Alternative Names	Rabbit IgG (Fc specific, Heavy Chain)
Gene ID	
SwissProt ID	
Immunogen	

Background

Research Area

Image Data



Western blot analysis of Reduced Rabbit IgG, H+L, diluted at 1:10000.