
Product Name: Epsilon Tubulin Mouse Monoclonal Antibody**Catalog #: AMM86163**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application


Dilution Ratio	IHC 1:100-1:200
Molecular Weight	55kDa

Antigen Information


Gene Name	Epsilon Tubulin
Alternative Names	Epsilon Tubulin
Gene ID	51175;
SwissProt ID	Q9UJT0
Immunogen	Synthetic Peptide

Background**Research Area**

Image Data



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Colon Carcinoma Tissue using Epsilon Tubulin Mouse mAb diluted at 1:200.