
Product Name: Luciferase Mouse Monoclonal Antibody**Catalog #: AMM86164**

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	WB
Reactivity	Other
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide, 0.5% protective protein and 50% glycerol.
Purification	Affinity Purification

Application


Dilution Ratio	WB 1:1000-1:2000
Molecular Weight	60kDa

Antigen Information

Gene Name	Luciferase
Alternative Names	Luciferase
Gene ID	
SwissProt ID	P08659
Immunogen	Synthetic Peptide

Background**Research Area**

Image Data



Western blot analysis of Firefly Luciferase Transfected Hela Cell Lysate using Luciferase Mouse mAb diluted at 1:2000.