

Product Name: Acetyl Lysine Mouse Monoclonal Antibody
Catalog #: AMM86166

Summary

Production Name	Acetyl Lysine Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Other

Performance

Conjugation	Unconjugated
Modification	Acetylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Acetyl Lysine
Alternative Names	Acetyl Lysine
Gene ID	
SwissProt ID	P68431.Purified Protein

Application

Dilution Ratio	WB:1:1000-1:2000,IHC:1:200-500
Molecular Weight	Calculated MW: 15 kDa; Observed MW: 15 kDa


Background

Product Name: Acetyl Lysine Mouse Monoclonal Antibody
Catalog #: AMM86166




Research Area

Image Data



Western blot analysis of 1) Mouse Brain Tissue, 2) Hela, 3) Hela+TSA Treated using Acetyl Lysine Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Acetyl Lysine Monoclonal Antibody.

Note

For research use only.