

**Product Name: Pan Methylated Lysine Mouse
Monoclonal Antibody
Catalog #: AMM86167**

Summary

Production Name	Pan Methylated Lysine Mouse Monoclonal Antibody
Description	Mouse Monoclonal Antibody
Host	Mouse
Application	WB,IHC
Reactivity	Other

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide and 50% glycerol.
Purification	Affinity Purification

Immunogen

Gene Name	Pan Methylated Lysine
Alternative Names	Pan Methylated Lysine
Gene ID	
SwissProt ID	P68431.Conjugated Protein

Application

Dilution Ratio	WB:1:1000-1:2000,IHC:1:200-500
Molecular Weight	Calculated MW: 15 kDa; Observed MW: 15 kDa


Background

**Product Name: Pan Methylated Lysine Mouse
Monoclonal Antibody
Catalog #: AMM86167**




Research Area

Image Data



Western blot analysis of Hela using Pan Methylated Lysine Monoclonal Antibody.



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma using Pan Methylated Lysine Monoclonal Antibody.

Note

For research use only.