
Product Name: Acetyl-NF kB P65(K314/K315) (2A11) Mouse Monoclonal Antibody
Catalog #: AMM86170

For research use only.

Summary

Description	Mouse monoclonal Antibody
Host	Mouse
Application	IHC
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Acetylated
Isotype	IgG
Clonality	Monoclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Purified antibody in PBS with 0.05% sodium azide,0.5%protective protein and 50% glycerol.
Purification	Affinity Purification

Application

Dilution Ratio	IHC 1:50-1:200
Molecular Weight	65kDa

Antigen Information


Gene Name	Acetyl NF kB P65(K314/K315) (2A11)
Alternative Names	Transcription factor p65 (Nuclear factor NF-kappa-B p65 subunit) (Nuclear factor of kappa light polypeptide gene enhancer in B-cells 3)
Gene ID	5970;
SwissProt ID	Q04206
Immunogen	Synthetic Peptide of Acetyl NF kB P65(K314/K315)

Background


Research Area

PI3K-Akt signaling pathway, MAPK signaling pathway

Image Data



Immunohistochemical analysis of paraffin-embedded Human Breast Carcinoma Tissue using Acetyl NF κB P65(K314/K315) Mouse mAb diluted at 1:200.



Immunohistochemical analysis of paraffin-embedded Human Lung Carcinoma Tissue using Acetyl NF κB P65(K314/K315) Mouse mAb diluted at 1:200.