

Summary

Production Name	Cbfa2t2 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Mouse, Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2c
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Cbfa2t2
Alternative Names	Cbfa2t2
Gene ID	
SwissProt ID	O70374. MTGR1 peptide

Application

Dilution Ratio	IHC:1:50
Molecular Weight	/


Background

Product Name: Cbfa2t2 Rat Monoclonal Antibody
Catalog #: AMR83163



Research Area

Image Data



IHC-P analysis of testis tissue by anti-Cbfa2t2 antibody (AMR83163)

Note

For research use only.