
Summary

Production Name	Chtop Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Chtop
Alternative Names	Chtop
Gene ID	
SwissProt ID	Q9CY57. The N-terminal fragment of mouse CHTOP

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: Chtop Rat Monoclonal Antibody
Catalog #: AMR83169



Research Area

Image Data



IHC-P analysis of mouse spleen tissue by anti- Chtop antibody (AMR83169)

Note

For research use only.