

Product Name: Lgi4 Rat Monoclonal Antibody
Catalog #: AMR83170

Summary

Production Name	Lgi4 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Lgi4
Alternative Names	Lgi4
Gene ID	
SwissProt ID	Q8K1S1. Purified recombinant fragment of mouse Lgi4

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: Lgi4 Rat Monoclonal Antibody
Catalog #: AMR83170



Research Area

Image Data



IHC-P analysis of mouse cerebellum tissue by anti- Lgi4 antibody (AMR83170)

Note

For research use only.