

Product Name: Ctf Rat Monoclonal Antibody
Catalog #: AMR83247



Summary

Production Name	Ctf Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB,IP
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Ctf
Alternative Names	Ctf
Gene ID	
SwissProt ID	Q61164.Purified recombinant fragment of mouse Ctf

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200
Molecular Weight	135kDa


Background

Product Name: Ctfc Rat Monoclonal Antibody
Catalog #: AMR83247




Research Area


Image Data



Lane 1: 15 μ g of HeLa lysate; Lane 2: 15 μ g of Jurkat lysate; Lane 3: 15 μ g of MOLT-4 lysate; Result: AMR83247 can detect CTCF by Western blotting;



Western blot analysis using CTCF mouse mAb against NIH/3T3 lysate



Lane 1: MOLT-4 lysate; Lane 2: CTCF immunoprecipitated from MOLT-4 lysate by AMR83247; Lane 3: The same as Lane 2 but use IgG isotype control antibody; Result: AMR83247 can immunoprecipitate CTCF;

Product Name: Ctf Rat Monoclonal Antibody
Catalog #: AMR83247



Note

For research use only.