

**Product Name: Ctcfl Rat Monoclonal Antibody**  
**Catalog #: AMR83249**

---

## Summary

<b>Production Name</b>	Ctcfl Rat Monoclonal Antibody
<b>Description</b>	Rat Monoclonal antibody
<b>Host</b>	Rat
<b>Application</b>	WB
<b>Reactivity</b>	Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Rat IgG2a
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	Ctcfl
<b>Alternative Names</b>	Ctcfl
<b>Gene ID</b>	
<b>SwissProt ID</b>	A2APF3. Purified recombinant fragment of mouse Ctcfl

## Application

<b>Dilution Ratio</b>	WB:1:500-1:1000
<b>Molecular Weight</b>	72kDa

## Background


**Product Name: Ctcfl Rat Monoclonal Antibody**  
**Catalog #: AMR83249**

---



## Research Area

## Image Data



Western blot analysis using mouse CTCFL mAb against mouse testis tissue lysate

## Note

For research use only.