

**Product Name: Smc3 Rat Monoclonal Antibody**  
**Catalog #: AMR83258**



## Summary

<b>Production Name</b>	Smc3 Rat Monoclonal Antibody
<b>Description</b>	Rat Monoclonal antibody
<b>Host</b>	Rat
<b>Application</b>	WB,IHC
<b>Reactivity</b>	Mouse

## Performance

<b>Conjugation</b>	Unconjugated
<b>Modification</b>	Unmodified
<b>Isotype</b>	Rat IgG2a
<b>Clonality</b>	Monoclonal
<b>Form</b>	Liquid
<b>Storage</b>	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
<b>Buffer</b>	Purified antibody in PBS with 0.05% sodium azide
<b>Purification</b>	Affinity Purification

## Immunogen

<b>Gene Name</b>	Smc3
<b>Alternative Names</b>	Smc3
<b>Gene ID</b>	
<b>SwissProt ID</b>	Q9CW03.Purified recombinant fragment of mouse Smc3

## Application

<b>Dilution Ratio</b>	WB:1:500-1:1000,IHC:1:100-1:200
<b>Molecular Weight</b>	55kDa


## Background

**Product Name: Smc3 Rat Monoclonal Antibody**  
**Catalog #: AMR83258**




## Research Area

## Image Data



Western blot analysis using SMC3 mAb against mouse testicles lysate and NIH/3T3 lysate



IHC-P analysis of colon tissue by anti-SMC3 antibody (AMR83258)

## Note

For research use only.