

Product Name: Stag2 Rat Monoclonal Antibody
Catalog #: AMR83262



Summary

Production Name	Stag2 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgM
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Stag2
Alternative Names	Stag2
Gene ID	
SwissProt ID	O35638. Purified recombinant fragment of mouse Stag2

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	130kDa


Background

Product Name: Stag2 Rat Monoclonal Antibody
Catalog #: AMR83262



Research Area

Image Data



Lane 1: 15 μ g of K-562 lysate; Lane 2: 15 μ g of Jurkat lysate; Lane 3: 15 μ g of MCF7 lysate; Result: AMR83262 can detect Stag2 by Western blotting;

Note

For research use only.