

Product Name: MAPRE3 Rat Monoclonal Antibody
Catalog #: AMR83268



Summary

Production Name	MAPRE3 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	MAPRE3
Alternative Names	MAPRE3
Gene ID	
SwissProt ID	Q9UPY8. Purified recombinant fragment of human MAPRE3

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	32kDa


Background

Product Name: MAPRE3 Rat Monoclonal Antibody
Catalog #: AMR83268



Research Area

Image Data



Western blot analysis using MAPRE3 mAb against mouse brain lysate

Note

For research use only.