

Product Name: Rab11b Rat Monoclonal Antibody
Catalog #: AMR83354



Summary

Production Name	Rab11b Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Rab11b
Alternative Names	Rab11b
Gene ID	
SwissProt ID	P46638. Purified recombinant fragment of mouse Rab11b

Application

Dilution Ratio	IHC:1:50
Molecular Weight	/


Background

Product Name: Rab11b Rat Monoclonal Antibody
Catalog #: AMR83354



Research Area

Image Data



IHC-P analysis of mouse brain tissue by Rab11b antibody (AMR83354)

Note

For research use only.