

Product Name: Histone H3K9ac Rat Monoclonal Antibody
Catalog #: AMR83415

Summary

Production Name	Histone H3K9ac Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB,IP,IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG1
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Histone H3K9ac
Alternative Names	Histone H3K9ac
Gene ID	
SwissProt ID	P68431/P84243. Purified recombinant fragment of human Histone H3K9ac

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200,IHC:1:500
Molecular Weight	17kDa


Background

Product Name: Histone H3K9ac Rat Monoclonal Antibody
Catalog #: AMR83415



Research Area

Image Data



Microtiter wells were coated with various peptides. AMR83415 was used as the primary antibody and peroxidase conjugated goat anti-rat IgG was used as the the secondary antibody.

Note

For research use only.