

Product Name: Cmtm6 Rat Monoclonal Antibody
Catalog #: AMR83416



Summary

Production Name	Cmtm6 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Mouse, Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Cmtm6
Alternative Names	Cmtm6
Gene ID	
SwissProt ID	Q9CZ69. Purified recombinant fragment of mouse Cmtm6

Application

Dilution Ratio	IHC:1:20-1:100
Molecular Weight	/


Background

Product Name: Cmtm6 Rat Monoclonal Antibody
Catalog #: AMR83416



Research Area

Image Data



IHC-P analysis of endometrium tissue by Cmtm6 antibody (AMR83416)

Note

For research use only.