

Product Name: TMEM161B Rat Monoclonal Antibody
Catalog #: AMR83431



Summary

Production Name	TMEM161B Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	TMEM161B
Alternative Names	TMEM161B
Gene ID	
SwissProt ID	Q8NDZ6. Purified recombinant fragment of human TMEM161B

Application

Dilution Ratio	IHC:1:100-1:200
Molecular Weight	/


Background

Product Name: TMEM161B Rat Monoclonal Antibody
Catalog #: AMR83431



Research Area

Image Data



IHC-P analysis of kidney tissue by TMEM161B antibody (AMR83431)

Note

For research use only.