

Product Name: MAPRE3 Rat Monoclonal Antibody
Catalog #: AMR83580



Summary

Production Name	MAPRE3 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	IP,IHC
Reactivity	Human, Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	MAPRE3
Alternative Names	MAPRE3
Gene ID	
SwissProt ID	Q9UPY8.Purified recombinant fragment of human MAPRE3

Application

Dilution Ratio	IP:1:100-1:200,IHC:1:500
Molecular Weight	/


Background

Product Name: MAPRE3 Rat Monoclonal Antibody
Catalog #: AMR83580




Research Area

Image Data



Lane 1: Mouse brain lysate ; Lane 2: MAPRE3 immunoprecipitated from mouse brain lysate by AMR83580 ; Lane 3: The same Lane 2 but use IgG isotype control antibody ;Result: AMR83580 can immunoprecipitate MAPRE3;



IHC-P analysis of lung tissue by MAPRE3 antibody (AMR83580)

Note

For research use only.