

Product Name: Gata3 Rat Monoclonal Antibody
Catalog #: AMR83582



Summary

Production Name	Gata3 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB,IHC
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Gata3
Alternative Names	Gata3
Gene ID	
SwissProt ID	P23772. Purified recombinant fragment of mouse Gata3

Application

Dilution Ratio	WB:1:500-1:1000,IHC:1:500
Molecular Weight	55kDa


Background

Product Name: Gata3 Rat Monoclonal Antibody
Catalog #: AMR83582



Research Area

Image Data



Lane 1: 15 µg of Jurkat lysate Lane 2: 15 µg of SH-SY5Y lysate Lane 3: 15 µg of MCF-7 lysate

Note

For research use only.