

Product Name: NIPBL Rat Monoclonal Antibody
Catalog #: AMR83586



Summary

Production Name	NIPBL Rat Monoclonal Antibody
Description	Rat monoclonal antibody
Host	Rat
Application	IHC
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	NIPBL
Alternative Names	NIPBL
Gene ID	
SwissProt ID	Q6KC79-2.Purified recombinant fragment of human NIPBL

Application

Dilution Ratio	IHC 1:200-1:500
Molecular Weight	/


Background

Product Name: NIPBL Rat Monoclonal Antibody
Catalog #: AMR83586



Research Area

Image Data



IHC-P analysis of stomach tissue by NIPBL antibody (AMR83586)

Note

For research use only.