

Product Name: Clasp1 Rat Monoclonal Antibody
Catalog #: AMR83591



Summary

Production Name	Clasp1 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB,IP
Reactivity	Mouse, Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2a
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Clasp1
Alternative Names	Clasp1
Gene ID	
SwissProt ID	Q80TV8.Purified recombinant fragment of mouse Clasp1

Application

Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200
Molecular Weight	180kDa


Background

Product Name: Clasp1 Rat Monoclonal Antibody
Catalog #: AMR83591




Research Area

Image Data



Lane 1: 15 μ g of U-937 lysate; Lane 2: 15 μ g of K-562 lysate; Lane 3: 15 μ g of rat heart tissue lysate; Result: AMR83591 can detect Clasp1 by Western blotting;



Lane 1: K-562 lysate ; Lane 2: Clasp1 immunoprecipitated from K-562 lysate by AMR83591 ; Lane 3: The same as Lane 2 but use IgG isotype control antibody ; Result: AMR83591 can immunoprecipitate Clasp1;

Note

For research use only.