

Summary

Production Name	Clasp2 Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB,IP
Reactivity	Mouse, Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Clasp2
Alternative Names	Clasp2
Gene ID	
SwissProt ID	Q8BRT1. Purified recombinant fragment of mouse Clasp2

Application


Dilution Ratio	WB:1:500-1:1000,IP:1:100-1:200
Molecular Weight	180kDa

Background

Product Name: Clasp2 Rat Monoclonal Antibody
Catalog #: AMR83599

Research Area

Image Data



Lane 1: 15 µg of HEK-293 lysate Lane 2: 15 µg of Jurkat lysate Lane 3: 15 µg of HeLa lysate

Note

For research use only.