

Product Name: Rab6b Rat Monoclonal Antibody
Catalog #: AMR83600



Summary

Production Name	Rab6b Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB
Reactivity	Mouse, Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG2b
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Rab6b
Alternative Names	Rab6b
Gene ID	
SwissProt ID	P61294. Purified recombinant fragment of mouse Rab6b

Application

Dilution Ratio	WB:1:500-1:1000
Molecular Weight	180kDa


Background

Product Name: Rab6b Rat Monoclonal Antibody
Catalog #: AMR83600



Research Area

Image Data



Lane 1: 15 μ g of SH-SY5Y lysate; Lane 2: 15 μ g of mouse brain tissue lysate; Lane 3: 15 μ g of rat brain tissue lysate; Result:
AMR83600 can detect Rab6b by Western blotting;

Note

For research use only.