

Product Name: Tmem176b Rat Monoclonal Antibody
Catalog #: AMR83673



Summary

Production Name	Tmem176b Rat Monoclonal Antibody
Description	Rat Monoclonal antibody
Host	Rat
Application	WB
Reactivity	Mouse

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	Rat IgG
Clonality	Monoclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Purified antibody in PBS with 0.05% sodium azide
Purification	Affinity Purification

Immunogen

Gene Name	Tmem176b
Alternative Names	Tmem176b
Gene ID	
SwissProt ID	Q9R1Q6. Purified recombinant fragment of mouse Tmem176b

Application

Dilution Ratio	WB:1:2000
Molecular Weight	26kDa,28kDa


Background

Product Name: Tmem176b Rat Monoclonal Antibody
Catalog #: AMR83673



Research Area

Image Data



Western blot analysis using Tmem176b mAb against mouse liver tissue lysate and mouse intestine tissue lysate

Note

For research use only.