

PBS Buffer (10×)

Catalog No.: PA10008

Basic Information

Product name	PBS Buffer (10×)
Sizes	500 mL
Storage	RT
Shipping	RT
Validity	24 months

Product Introduction

PBS Buffer (10×) is a 10-fold concentrated phosphate-buffered saline widely used as a balanced salt solution in biological and biochemical research. After 10-fold dilution it becomes 1× PBS and can be used for IHC, IF, WB, and other applications.

This product has been sterile-filtered through 0.1 µm.

Working Solution Preparation

Add 100 mL of the 10× concentrate to 900 mL distilled water and mix well to obtain 1× PBS containing 0.01 M phosphate at pH 7.3–7.5.

Notes

1. This product is a 10×d concentrate; dilute as required before use.
2. Wear a lab coat and disposable gloves for your safety and health.
3. For research use only; not intended for clinical diagnosis or therapy.

This product is for research use only!