

PBST Buffer (10×)

Catalog No.: PA10009

Basic Information

Product name	PBST Buffer (10×)
Sizes	500 mL
Storage	RT
Shipping	RT
Validity	12 months

Product Introduction

This 10× PBST buffer is specially formulated for high-sensitivity assays such as IHC, IF, and WB. Compared with PBS, it removes non-specifically bound antibodies more efficiently and significantly lowers background noise.

Working Solution Preparation

Add 10 mL of the 10× concentrate to 90 mL distilled water and mix well.

Notes

1. This product is not manufactured under sterile conditions.
2. Use promptly after opening; alternatively, store at 4 °C.
3. This product is potassium-free.
4. Wear a lab coat and disposable gloves for your safety and health.
5. For research use only; not intended for clinical diagnosis or therapy.

This product is for research use only!