

TBS Buffer (10×)

Catalog No.: PA10010

Basic Information

Product name	TBS Buffer (10×)
Sizes	500 mL
Storage	RT
Shipping	RT
Validity	12 months

Product Introduction

Tris-buffered saline (TBS) is used for washing and for preparing blocking solutions during immunohistochemistry (IHC), immunofluorescence (IF), in-situ hybridization (ISH), ELISA, Western blot (WB), and other immunoassays.

This product is a 10× concentrate and must be diluted 10-fold before use.

Working Solution Preparation

Add 100 mL of the 10× concentrate to 900 mL distilled water and mix well to obtain 1× TBS working solution.

Notes

1. This product is not manufactured under sterile conditions.
2. After dilution to 1× working solution, store at 2–8 °C and use within one week.
3. Wear a lab coat and disposable gloves for your safety and health.
4. For research use only; not intended for clinical diagnosis or therapy.

This product is for research use only!