

Product Name: PRC1 Rabbit Polyclonal Antibody**Catalog #: APRab00002**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,FC,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,FC 1:50-1:100,IP 1:20-1:50
Molecular Weight	Calculated MW: 72 kDa; Observed MW: 72 kDa

Antigen Information

Gene Name	PRC1
Alternative Names	Protein regulator of cytokinesis 1
Gene ID	9055
SwissProt ID	O43663
Immunogen	A synthetic peptide corresponding to target protein


Background

Cross-links antiparrallel microtubules at an average distance of 35 nM. Essential for controlling the spatiotemporal formation of the midzone and successful cytokinesis. Required for KIF14 localization to the central spindle and midbody.

Research Area

Cell Biology

Image Data



Western blot analysis of PRC1 in HeLa lysates using PRC1 antibody.