

Product Name: Thrombomodulin Rabbit Polyclonal Antibody
Catalog #: APRab00023

Summary

Production Name	Thrombomodulin Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,IHC-P,ICC/IF,IP
Reactivity	Human

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	THBD
Alternative Names	CD141; Fetomodulin; THBD; THRM; thrombomodulin; TM
Gene ID	7056
SwissProt ID	P07204.

Application

Dilution Ratio	WB: 1:500-1:1000 IHC: 1:50-1:100 IF: 1:50-1:200 IP: 1:20
Molecular Weight	Calculated MW: 60 kDa; Observed MW: 60 kDa

Product Name: Thrombomodulin Rabbit Polyclonal Antibody
Catalog #: APRab00023




Background

Thrombomodulin (TM), also called CD141, is a type I membrane receptor that is specific to endothelial cells. TM has a cysteine-rich extracellular domain with six EGF-like regions. It forms a complex with Thrombin, which activates Protein C to generate activated Protein C (APC), an anticoagulant enzyme. APC together with Protein S inhibits coagulation by inactivating Factors Va and VIIIa. Deletion of the TM gene results in embryonic lethality in mice.


Research Area

Cardiovascular

Image Data



Western blot analysis of Thrombomodulin in Human placenta lysates using Thrombomodulin antibody.



Immunohistochemistry analysis of paraffin-embedded Human bladder using Thrombomodulin antibody. High-pressure and temperature Sodium Citrate pH 6.0 was used for antigen retrieval.

Note

For research use only.