

Product Name: Cyclophilin B Rabbit Polyclonal Antibody**Catalog #: APRab00033**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000
Molecular Weight	Calculated MW: 24 kDa; Observed MW: 18 kDa

Antigen Information

Gene Name	PPIB
Alternative Names	PPIB; CYPB; Peptidyl-prolyl cis-trans isomerase B; PPIase B; CYP-S1; Cyclophilin B; Rotamase B; S-cyclophilin; SCYLP
Gene ID	5479
SwissProt ID	P23284
Immunogen	A synthetic peptide of human Cyclophilin B

Background


PPIases accelerate the folding of proteins. It catalyzes the cis-trans isomerization of proline imidic peptide bonds in

oligopeptides.

Research Area

Immunology

Image Data



Western blot analysis of Cyclophilin B in HepG2 lysates using Cyclophilin B antibody.