

Product Name: Rab9 Rabbit Polyclonal Antibody
Catalog #: APRab00049



Summary

Production Name	Rab9 Rabbit Polyclonal Antibody
Description	Rabbit Polyclonal Antibody
Host	Rabbit
Application	WB,ICC/IF
Reactivity	Human,Mouse,Rat

Performance

Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Storage	Store at 4°C short term. Aliquot and store at -20°C long term. Avoid freeze/thaw cycles.
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Immunogen

Gene Name	RAB9A
Alternative Names	RAB9A; DmRab9; Ras related protein Rab 9A
Gene ID	9367
SwissProt ID	P51151.

Application

Dilution Ratio	WB: 1:500-1:1000 IF: 1:50-1:200
Molecular Weight	Calculated MW: 23 kDa; Observed MW: 23 kDa

Background

Involved in the transport of proteins between the endosomes and the trans Golgi network. Interacts (preferentially in its

Product Name: Rab9 Rabbit Polyclonal Antibody
Catalog #: APRab00049




GTP-bound form) with GCC2 (via its GRIP domain).

Research Area

Tags & Cell Markers

Image Data



Western blot analysis of Rab9 in HepG2 lysates using Rab9 antibody.

Note

For research use only.