
Product Name: TFG Rabbit Polyclonal Antibody**Catalog #: APRab00119**

For research use only.

Summary

Description	Rabbit polyclonal Antibody
Host	Rabbit
Application	WB,IHC,ICC/IF,IP
Reactivity	Human,Mouse,Rat
Conjugation	Unconjugated
Modification	Unmodified
Isotype	IgG
Clonality	Polyclonal
Form	Liquid
Concentration	1mg/ml
Storage	Aliquot and store at -20°C (valid for 12 months). Avoid freeze/thaw cycles.
Shipping	Ice bags
Buffer	Rabbit IgG in phosphate buffered saline , pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
Purification	Affinity Chromatography

Application

Dilution Ratio	WB 1:500-1:1000,IHC 1:50-1:100,ICC/IF 1:50-1:200,IP 1:20-1:50
Molecular Weight	Calculated MW: 43 kDa; Observed MW: 55 kDa

Antigen Information

Gene Name	TFG
Alternative Names	TF6; HMSNP; SPG57; TRKT3
Gene ID	10342
SwissProt ID	Q92734
Immunogen	A synthetic peptide of human TRK fused gene


Background

Plays a role in the normal dynamic function of the endoplasmic reticulum (ER) and its associated microtubules.

Research Area

Immunology

Image Data



Western blot analysis of TRK fused gene in HeLa lysates using TFG antibody.